



COMPLETE PANEL OF TESTS

ADAPTED TO CHILDREN

The **VT1 Junior Edition** offers a large choice of tests to screen the main visual disorders in under 5 minutes. These tests are specifically designed for children > 5 years old.

Tests available

Acuity

(20/400 20/200, 20/100 20/50 20/40 20/30 20/25 20/20) Symbols for children > 5 years old SLOAN letters for children \geq 7 years old Hyperopia +1 δ

Low vision with symbols and SLOAN letters

- Color test
- Vertical & Lateral Phorias
- **Depth Perception** (800" 400" 200" 100" 80" 40")
- Fusion
- · Visual Field ((horizontal, nasal & vertical)

1 SROZC ODVNZ NZHOD AND Test Examples

ALL INCLUSIVE PACKAGE

WITH LCD REMOTE AND SOFTWARE

- ✓ Drive your device from your computer or remote control at your convenience
- ✓ Use the LCD remote to perform tests when your laptop is unavailable. (Response forms provided with the CD)
- ✓ Use the software to benefit from all the advantages of a computerized device (fast, results export, interfaced with EMRs, sequences customization, etc.)

CREATE YOUR OWN

TEST SEQUENCES



To save time when performing tests, customize your own test sequences using either the software or the remote control (on request).

2 test sequences are already included for your convenience.



DEPISTEOInnovation

DRIVE YOUR
DEVICE WITH
EASE!







WHEELED CARRYING CASE

- ✓ Anti-shock wheels
- ✓ Reinforced interior
- ✓ Height-adjustable handle

ADDITIONAL ACCESSORIES



• Single-use face cover to add hygienic measure for the patient



 VT1 plastified dust cover protection



SUPPORT

Depisteo LLC developed guidelines to



YOUR DISTRIBUTOR

WARRANTY 2 years DIMENSIONS WEIGHT LIGHTING CONDITIONS

19,7 x 10,63 x 9,8" / 50x27x25cm 4.8 kg / 10.7 lbs
• High photopic (160 candelas) (NOMINAL VALUES) Low photopic (80 candelas) Lit by 16 white LEDs driven by microcontroller and diffuser system LIGHTING SYSTEM

DISTANCES • Far Vision: (5.0 ± 0.1) m / (20.0 ± 0.4) ft

 Intermediate Vision: (60.0 ± 0.5) cm/ (24.0 ± 0.2)" Near Vision: (35.5 ± 0.25) cm/ (14.00 ± 0.1)"

 Hyperopia Lens: +I dioptre
High definition photographic film TEST BAND REACTIVITY

Time between two neighbouring tests: 700ms
Time to go from one slide to another: 1s **TEST TIME** · Routine test: 3 minutes In-depth test: 5 minutes

OPERATING TEMPERATURE between 15 and 35°C / between 59 and 95°F HUMIDITY 75% maximum POWER SUPPLY 12 VDC from an external medical power supply POWER ABSORBED

REFERENCE STANDARDS EN 60601-1, EN 60601-1-2, IEC 60601-1-6, EN 62366-1, EN ISO 10993-5, EN ISO 10993-10, NF EN ISO 14971, NF EN

62304/AI NF EN 1041+AI, EN ISO 15223-1, ISO 8596, ANSI Z80.21, NF EN ISO 15004-2 :2007 § 5.4.1.6, EN ISO 10993-1

MARKING MEDICAL CLASS Class I APPLIED PART Type B



550 Pharr Road, NE, Suite 201 Atlanta GA 30305

Phone: +1 800 496 8290 Fax: + 1 800 496 8290 contact@depisteo.com

 $C \in$