

VISUAL SCREENING



VISIOCLICK[®]

AUTOMATION AT YOUR FINGERTIPS

BENEFITS

- Adapts to your VT1[®] fleet
- Automation of all tests
- Multilingual voice instructions



CE



Convert your VT1® Master or Master GT to an automated device. VisioClick also adapts to previous versions (from 2012).

Compact, lightweight, and operational from the first connection. A simple update of the VT1® software is all it takes.

No operator presence required during tests. Multilingual voice instruction ensures global understanding.

HOW DOES IT WORK ?

4-step tests. No operator intervention

STEP 1: HEADPHONES

The patient puts on the headphones and receives test instructions displayed in the VT1®. Multilingual Voice Instructions available.

STEP 2: RESPONSE UNIT

Patient responds directly to instructions by clicking on his unit.

STEP 3: RECORDING

Results are automatically saved.

STEP 4: END OF TEST

Information is given to the patient and operator indicating the test is completed. Results are recorded and can be interfaced with all «business software».

NEW !



EARPAD COVERS PACK

Single-use biocompatible for all your tests.

VisioClick® delivered with a pack of earpad covers.

Benefits:

- Increases lifespan of earpads
- No sound reduction
- A more hygienic solution



VT1® SOFTWARE

- Windows XP SP3, Windows Vista, Windows 7, Windows 8 and Windows 10 compatible.
- Interfacing with market software packages (OHM, Cority...)

Specific software features:

- Device can perform all tests individually in the order of your choice or a pre-programmed entire test protocol
- Creation of customizable profession-related sequences that can be performed in automatic mode
- Multilingual voice instructions (French, English, German, Spanish, Dutch, Italian, Polish, Portuguese, Turkish, Arabic, etc.)

YOUR DISTRIBUTOR

WARRANTY	2 years (without headphones support)
SIZE	9x5x2 inch
WEIGHT	1.55 lbs
CONNECTIVES	Jack 0.12 inch (to audio headphones), USB Hub (to VT1® device)
POWER SUPPLY VOLTAGE	5VDC (via USB)
REFERENCE STANDARDS	EN60601-1, EN60601-1-2
MARKING	CE
MEDICAL CLASS	Classe I (directive 93/42/CEE)



550 Pharr Road, Suite 201
Atlanta GA 30305

Phone: +1 800 496 8290
Fax: +1 800 496 8290
contact@depisteo.com

